


Corneal Transplant



100% sterile solution*

The Moria system for creating ultra-thin and conventional corneal lamellar grafts

*hand-piece and all accessories either autoclavable or disposable

Moria
Leading Innovations in Ophthalmology

CB Single-Use Heads




- Range of heads:
 - 350-micron
 - 300-micron
 - 250-micron
 - 200-micron
 - 130-micron
 - 110-micron
 - 90-micron
 - 50-micron
- Based on corneal thickness, can be used in various combinations to create ultra-thin lamellar grafts.
- Can also be used individually to create conventional thickness grafts.

Single-Use Artificial Chamber




- Rotating cover allows initiation of cuts in opposite quadrants of the cornea, enabling creation of grafts that are more planar than conventional grafts.
- Creates the largest possible diameter graft.
- Can also be customized for other procedures in which a microkeratome is not used.

Inflation Tubing Set



- Consistently maintains pressure required for creation of accurate, predictable ultra-thin grafts.
- Enables use of two different methods:
 - with corneal storage medium, or
 - with BSS and cohesive viscoelastic.


CB Micro-keratome



- Rotating CB turbine
- Evolution 3E control unit

Worldwide, most corneal lamellar grafts are created with the Moria system.

The global market leader in equipment for creating corneal lamellar grafts

Product number	Description		Packaging
19182	Single-Use Artificial Chamber (alone)		1
19178/...	CBSU heads (50-90-110-130-200-250-300-350)		5
19180	CBSU guide ring and blue screw		5
19181	Infusion tubing set for Artificial Chamber		10

To obtain more information



MORIA SA

15, rue Georges Besse
92160 Antony
FRANCE
Phone: +33 (0) 1 46 74 46 74
Fax: +33 (0) 1 46 74 46 70
www.moria-surgical.com

MORIA Inc

1050 Cross Keys Drive
Doylestown, PA 18902
USA
Phone: (800) 441 1314
Fax: +1 (215) 230 7670
www.moria-surgical.com

MORIA JAPAN KK

Inami Bldg. 6F
3-24-2 Hongo Bunkyo-ku
Tokyo 113-0033
JAPAN
Phone: +81 3 5842 6761
Fax: +81 3 5842 6762
www.moriajapan.com

MORIA in CHINA

Moria Commercial (China) Co. Ltd
Rm H, I, 6FL, Kai Li Building
N° 432 Huai Hai West Road
Shanghai 200052
CHINA
Phone: +86 21 5258 5066
Fax: +86 21 5258 5067
www.moria-surgical.com.cn

Moria
Leading Innovations in Ophthalmology